



## Siemens Motion™

The hearing instrument for people  
who are too busy living.

[www.usa.siemens.com/hearing](http://www.usa.siemens.com/hearing)

**SIEMENS**



Motion is ergonomically designed for a comfortable fit.

Whether you're working in the garden or spending quality time with friends, you're too busy to let a hearing loss get in your way. That's why we developed Siemens Motion™. It's an easy, comfortable solution with all the technology you need to stay connected to your world. Siemens Motion. Because you have better things to do.



## Easy does it.

Siemens created Motion for hassle-free hearing. It's fully automatic, so there's no need to change programs or adjust volume. Motion 700 even has a SoundLearning™ feature that remembers your volume, bass and treble preferences for a consistent and more natural hearing experience. Designed for easy handling and simple operation, Motion is available with a choice of easy-to-use optional remote controls. It also features Autophone® technology for seamless performance when using the phone. Motion BTEs are rechargeable, so you don't have to fiddle with batteries. For greater convenience and versatility, they also work with regular hearing instrument batteries. What could be easier?



To recharge your Motion BTE instruments, simply place them in the easy-to-use charger — up to two instruments at a time. After just five hours of charging, the instruments are ready to go all day and evening.

## Act naturally.

Knowing where sounds are coming from is critical to a natural hearing experience. To help you accurately locate the source of sounds, Motion 700 BTEs combine two innovative technologies. TruEar® simulates the function of your outer ear to help you distinguish sounds in front of you from those behind you. And e2e wireless® 2.0 synchronizes your instruments, so they work together to help you identify the source of sounds from left to right. Motion 700 models also feature SoundBrilliance™; a technology designed to enhance your perception of higher frequencies for a much richer sound quality.



Motion is available with a choice of optional remote control devices, the ultra-slim ePen™ (left) and the ProPocket™ (right). Both provide fingertip control over virtually all functions.

## All the right connections.

If you enjoy watching TV or listening to music, you'll love Tek™, the optional enhancement that lets you wirelessly connect your Motion hearing instruments to your favorite audio devices. Tek works with all Motion 700 and 500 models through a dedicated remote— Tek™ Connect. Stylish and powerful, Tek Connect transforms your Motion hearing instruments into a high-performance stereo headset that lets you stream audio from your MP3 player, television, or other audio sources directly to your Motion instruments. That means you can listen to the TV at your preferred volume without disturbing others. Tek even connects you to your mobile phone with Bluetooth® wireless technology.

**TEK**) Siemens Motion is Tek-ready.

## Simply brilliant.

Motion offers a range of advanced technologies to provide unparalleled hearing performance.

### Speech and Noise Management

Motion automatically identifies and reduces background noise — like the buzz of traffic or the hum of a vacuum — without affecting important speech signals, so you can enjoy easier listening in any situation without excessive strain or fatigue.

### SoundSmoothing®

Unexpected noises like rustling paper and clanging dishes can be disturbing for hearing instruments wearers. That's why Motion is equipped with SoundSmoothing, which automatically softens sudden, annoying noises without affecting the quality of speech.

### FeedbackBlocker™

Motion is equipped with our advanced FeedbackBlocker. It searches for feedback patterns in sound frequencies and automatically reduces or eliminates annoying whistling.

### Automatic Adaptive Directional Microphone

Motion features an Automatic, Multichannel Adaptive Directional Microphone that's clinically proven to improve speech clarity in noisy situations like parties or restaurants where many sounds compete for your attention.


### eWindscreen®

On models with directional microphone, Motion's advanced eWindscreen feature monitors your environment for wind frequencies and then automatically softens them for improved comfort.


### Make your own motion

Designed for wearers with moderate through severe hearing losses, Motion is available in a wide range of BTE and custom models. A wide selection of colors is also available. To learn more, contact your Audiologist or Hearing Care Provider.

### Standard Colors



### Hair Tone Colors



### New Colors

# Your Audiologist or Hearing Care Provider:

The *Bluetooth*® word mark and logos are owned by the Bluetooth SIG, Inc., and any use of such marks by Siemens AG is under license. Other trademarks and trade names are those of their respective owners.

Hearing instruments help many people hear better, but no hearing instrument can solve every hearing problem nor restore normal hearing.

Information in this brochure is subject to change without notice. Copyright © 2008 Siemens Hearing Instruments, Inc. All rights reserved.

## Local Contact Information

Siemens Hearing Instruments, Inc.  
P.O. Box 1397  
Piscataway, NJ 08855-1397  
Phone: (800)766-4500  
Fax: (732)562-6696  
[www.usa.siemens.com/hearing](http://www.usa.siemens.com/hearing)

## Global Siemens Headquarters

Siemens AG  
Wittelsbacherplatz 2  
80333 Muenchen  
Germany

## Legal Manufacturer

Siemens Audiologische  
Technik GmbH  
Gebberstrasse 125  
91058 Erlangen  
Germany

## Global Business Unit

Siemens Audiologische  
Technik GmbH  
Gebberstrasse 125  
91058 Erlangen  
Germany  
Phone: +49 9131 308 0

## Global Siemens Healthcare Headquarters

Siemens AG  
Healthcare Sector  
Henkestrasse 127  
91058 Erlangen  
Germany  
Phone: +49 9131 84 0  
[www.siemens.com/healthcare](http://www.siemens.com/healthcare)



Better hearing. Better living.

Offices in:  
High Point, Greensboro, and Asheboro, NC  
<http://www.thehearingclinic.com/>  
[info@thehearingclinic.com](mailto:info@thehearingclinic.com)

**1-800-651-8551**

[www.siemens.com/hearing](http://www.siemens.com/hearing)